


CloseFit™

Narrow Implants



UNP CloseFit™
 Ø 2.75mmD
 Ultra Narrow Platform

NP CloseFit™
 Ø 3.0mmD
 Narrow Platform

* Optimally designed for narrow ridges, minimize the need of bone augmentation

* **OsseoFix™**

* **Double lead threads**
 2.4mm (2x1.2)
 For fast efficient insertion

* **Platform Switching**

* **Micro Threads:**
 Engineered to enhance bone regeneration, prevent bone loss & increase bone implant contact

Triple lead thread with 3.0mm (3x1)
 For fast efficient insertion

* **Length**

L-10mm | L-11.5mm | L-13 mm
 L-16mm | L-18 mm

*The features related to NP & UNP Implants

CloseFit™


Narrow Implants


With more than 10 years of proven clinical experience, these uniquely designed small diameter implants of 2.75mm & 3.0mm, have an especially narrow body & collar to provide a perfect solution that addresses narrow ridges in atrophic jaws for exceptional esthetic's results.

The strong Conical Connection designed to minimize micro movements.

Adin portfolio includes all components needed to support restorations on implant level, multi-unit level, cement retained and screw retained restorations, from single unit to full arch restoration, for temporary and permanent rehabilitation.


Our components gives a full CAD CAM solution for digital work flow, and are available in CAD software.





NP CloseFit™

Narrow Platform



Recommended insertion torque up to 50 Ncm


CloseFit™ NP	Ø	L	D3	Catalog No.
	3.0	10.0	2.3	ISPNP1030
		11.5	2.2	ISPNP1130
		13.0	2.2	ISPNP1330
		16.0	2.15	ISPNP1630
		18.0	2.1	ISPNP1830

Radiographic Template NP CloseFit™ | ASD0611

Drill Sequence

Recommended drilling speed up to 1200 rpm

Implant Diameter $\varnothing 3.0$ $\varnothing 3.0$



All measurements are in mm.


Cover Screw & Healing Abutments Recommended Torque 15(Ncm)

NP CloseFit™ Conical Cover Screw



NP0073

NP CloseFit™ Slim Healing Abutment



NP0046, GH2/L6.6


NP0047, GH3/L7.6

NP0048, GH4/L8.6

NP0049, GH5/L9.6

NP0050, GH6/L10.6

NP CloseFit™ Healing Abutment



NP0004, GH2/L6.5


NP0008, GH3/L7.5

NP0009, GH4/L8.5

NP0012, GH5/L9.5

NP0013, GH6/L10.5

NP CloseFit™ Conical Healing Abutment



NP0090, GH2/L6.9

NP0091, GH3/L7.9

NP0092, GH4/L8.9


NP0093, GH5/L9.9

NP0094, GH6/L10.9




NP™ Impression

Recommended Torque 15(Ncm)





Our components for taking impressions for all restoration types can be successfully utilized for digital or conventional protocols either implant-level or abutment level procedures.

Our broad range of solutions and flexible designs enable easy and accurate restorative planning, whether you're dealing with a single tooth, multiple teeth, bridge or full mouth rehabilitation.




Transfers		L	Ø	D1	H	GH	RH	Catalog No.
	NP Slim Closed Tray Transfer	15.2	3.2	2.5	13.0	-	-	NP5113
	NP Closed Tray Transfer	15.2	4.3	2.5	13.0	-	-	NP5213
	NP Slim Open Tray Transfer	18.0	5.0	2.5	14.0	-	-	NP5614
	NP Open Tray Transfer	18.0	6.0	2.5	14.0	-	-	NP5514
	NP Scan Body	9.2	-	2.5	7.0	-	-	NP0088
		12.2	-	2.5	10.0	-	-	NP0089

Analog/Digital Analog		L	Ø	D1	H	GH	RH	Catalog No.
	NP Implant Analog	12.6	3.0	2.5	-	-	-	NP0007
	NP digital analog 15(Ncm)	10.0	3.55	3.4	-	-	-	NP1000



Cement Retained Abutments

Recommended Torque 35(Ncm)





Abutments	L	Ø	D1	D2	H	GH	RH	Catalog No.
 NP Straight Slim Abutment No End Line	11.0	3.5	2.5	3.0	8.8	3.2	5.7	NP0010
 NP Straight Abutment	9.7	4.4	2.5	3.0	7.5	1.0	6.5	NP0019
	10.7	4.4	2.5	3.0	8.5	2.0	6.5	NP0011
	11.7	4.4	2.5	3.0	9.5	3.0	6.5	NP0064
 NP Angled Abutment 15° No End Line	11.7	4.3	2.5	1.9	9.5	1.4	8.2	NP0018
 NP Angled Abutment 15°	9.7	4.8	2.5	2.4	7.5	1.0	6.5	NP0015
	10.7	4.8	2.5	2.4	8.5	2.0	6.5	NP0026
	11.7	4.8	2.5	2.4	9.5	3.0	6.5	NP0063
 NP Angled Abutment 25° No End Line	10.0	3.8	2.5	1.9	7.8	1.6	6.5	NP0025
 NP Engaging Temporary Abutment	11.0	3.7	2.5	3.2	8.8	0.7	8.0	NP0043
 NP Non-Engaging Temporary Abutment	11.0	3.7	2.5	3.2	8.8	0.7	8.0	NP0044
 NP Engaging Plastic Cylinder Abutment (Casting)	11.0	3.7	2.5	3.0	8.8	0.7	8.0	NP0045
 NP Engaging Titanium / Plastic Abutment (Casting)	13.2	4.2	2.5	4.2	11.0	1.0	10.0	NP3095
	13.2	4.2	2.5	4.2	11.0	1.0	10.0	NP3097



TMA™ - Multi Unit

Recommended Torque 35(Ncm)









Our TMA™ (Trans Mucosal Abutment) system is indicated for multiple-unit, screw-retained restorations and may be used in combination with an implant level framework design.

The TMA™ system is used to elevate the restorative seating platform when restoration at implant level is not indicated or practical due to the depth or angle of the implant.

The TMA™ System is designed to allow better prosthetics access when using implants with divergent angles accommodate full and partial edentulous arches, especially when using angled implants and all-on-x technique for full arch restoration.

The TMA™ abutment system is available in straight or angled (17°, 30°, 45°), engaging the Conical Connection for indexing, with a selection of collar heights. Abutments are delivered with handles for easy insertion and to assist with seating.



TMA		L	Ø	D1	D2	H	GH	RH	Catalog No.
	NP Straight Trans Mucosal Abutment	4.2	4.9	2.5	-	2.0	1.0	2.10	NP0032
		5.2	4.9	2.5	-	3.0	2.0	2.10	NP0033
		6.2	4.9	2.5	-	4.0	3.0	2.10	NP0034
		7.2	4.9	2.5	-	5.0	4.0	2.10	NP3729
		8.2	4.9	2.5	-	6.0	5.0	2.10	NP3730
	NP Angled Trans Mucosal Abutment 17°	6.8	4.9	2.5	-	-	2.8	2.15	NP0035
	NP Angled Trans Mucosal Abutment 17°	7.9	4.9	2.5	-	-	3.9	2.15	NP0038
	NP Angled Trans Mucosal Abutment 30°	6.9	4.9	2.5	-	-	3.5	2.15	NP0036
	NP Angled Trans Mucosal Abutment 30°	8.1	4.9	2.5	-	-	4.9	2.15	NP0037

TMA components		L	∅	D1	D2	H	GH	RH	Catalog No.
	TMA Healing Cap	5.5	4.9	-	3.5	-	-	-	RS5005
		7	4.9	-	3.5	-	-	-	RS5003
	TMA Impression Coping Open Tray	11.0	4.9	-	3.1	-	-	-	RS5006
	TMA Closed Tray Transfer	8.0	4.9	-	4.0	-	-	-	RS5011
	TMA Scan Body	6	4.9	4.9	3.8	4.9	-	-	TM0006
		9	4.9	4.9	3.8	4.9	-	-	TM0009
	TMA Titanium Temporary Cylinder	12.0	4.9	-	3.5	-	1.5	10.5	RS4900
	TMA Cementing Cone	4.9	4.9	-	2.9	-	1.5	3.4	TM0010
	TMA Plastic Casting Sleeve w/ Prosthetic Screw	12.0	4.9	-	3.3	-	1.5	10.5	RS5001
	TMA Abutment Replica	16.8	4.9	-	1.7	-	-	-	RS5004
	TMA Digital Model Replica Recommended torque 15(Ncm)	12	4.9	-	3.8	-	-	-	RS1001




Overdenture Restoration

Recommended Torque 35(Ncm)



Ball Attachment		L	Ø	D1	D2	H	GH	RH	Catalog No.
	NP Ball Attachment	8.8	2.9	2.5	2.5	4.1	0.5	3.5	NP0020
		9.2	2.9	2.5	2.5	4.5	1.0	3.5	NP0021
		10.2	2.9	2.5	2.5	5.5	2.0	3.5	NP0022
		11.2	2.9	2.5	2.5	6.5	3.0	3.5	NP0023
		12.2	2.9	2.5	2.5	7.5	4.0	3.5	NP0024



Ball Attachment caps		L	Ø	D1	D2	H	GH	RH	Catalog No.
	Plastic Ball Cap White (Nylon)	2.60	4.10	-	3.60	-	-	-	RS2660
	Plastic Ball Cap Pink Soft (Nylon)	2.60	4.10	-	3.60	-	-	-	RS2662
	Titanium Ball Cap	2.90	4.70	-	4.70	-	-	-	RS2670
	Stainless Steel Ball Cap	2.90	4.70	-	4.70	-	-	-	RS2675SS

Ball Attachment		L	Ø	D1	D2	H	GH	RH	Catalog No.
	Ball Attachment Abutment Replica	16.60	4.00	-	2.5	-	0.80	2.7	RS5740



Screw - Retained

Recommended Torque 35(Ncm)



Flat Connection	L	Ø	D1	D2	H	GH	RH	Catalog No.
NP Flat Connection Abutment	7.3	4.5	2.5	-	2.0	2.0	-	NP0030
NP Flat Connection Abutment	8.8	4.5	2.5	-	2.0	4.0	-	NP0031
Flat Connection Temporary Cylinder Sleeve	11.5	4.5	-	3.9	11.1	0.4	10.7	FC5515
Flat Connection Implant Replica	12.6	4.5	-	3.4	-	-	-	FC5737
Flat Connection Castable Abutment for Centering Ring	13.6	4.6	-	3.0	13.6	0.1	13.5	FC6011
Flat Connection Titanium Retaining Screw	5.0	2.5	-	-	-	-	-	FC6015
Flat Connection Titanium Gluing Ring	3.8	4.5	-	3.1	3.4	0.6	2.8	FC6013
Flat Connection Long Laboratory Screw	17.0	2.5	-	-	-	-	-	FC6017
Flat Connection Centering Ring	3.8	4.5	-	3.2	3.4	0.6	2.8	FC6018
Flat Connection Transfer Impression Coping	11.0	4.5	-	4.2	-	-	10.2	FC5512

Screws		L	Ø	Catalog No.
	NP Abutment Screw	7.9	2.0	NP0001
	TMA Prosthetic Screw 15(Ncm)	4.1	2.1	RS3404
	NP Retaining Screw for Angled TMA	7.1	2.0	NP0042
	NP Retaining Screw for TMA 1mm	8.0	3.1	NP0039
	NP Retaining Screw for TMA 2mm	9.0	3.1	NP0040
	NP Retaining Screw for TMA 3mm	10.0	3.1	NP0041
	NP Retaining Screw for TMA 4mm	11.0	3.1	NP3409
	NP Retaining Screw for TMA 5mm	12.0	3.1	NP3412
	NP Open Tray Transfer Screw	28.0	3.1	NP5014
	NP Closed Tray Transfer Screw	20.0	2.8	NP5127
	NP Retrieval Screw for Abutments	13.4	2.0	NP0014



Hand & Torque Drivers

Recommended Torque 35(Ncm)



Star Hand Driver
23.4mm RP0064

*For Angeld TMA/TMA prosthetic screw



Star Hand Driver
23.4mm RP0004

*For Angeld TMA/TMA prosthetic screw



TMA/FC Female
Hand Hex Driver
22.4mm RS6193

*For straight TMA/Flat connection



TMA/FC Female
Hand Hex Driver
28.4mm RS6191

*For straight TMA/Flat connection



Wrench Ratchet
Surgical Instruments
Handle 16.7mm RA0050



UNP/NP Abutment
insertion tool 25.9mm UNP1001



NP 2mm Hex
Torque Driver
17mm NP0006



NP 2mm Hex
Torque Driver
27mm NP0005



Star Torque Driver
20.4mm RP0065

*For Angeld TMA/TMA prosthetic screw



Star Torque
Driver
25.4mm RP0005

*For Angeld TMA/TMA prosthetic screw



Torque Driver for
Ball Attachment
24.4mm RS6090



TMA/FC Hex Torque
Driver 20mm RS6197

*For straight TMA/Flat connection



TMA/FC Hex
Torque Driver
25mm RS6196

*For straight TMA/Flat connection



Handpiece Drivers

Recommended Torque 35(Ncm)



NP 2mm Hex Handpiece Adapter 23.7mm

NP0068



NP 2mm Hex Handpiece Adapter 27.7mm

NP0071



NP 2mm Hex Handpiece Adapter 36mm

NP0074



NP 2mm Hex Handpiece Adapter 20mm

NP0062



NP 2mm Hex Handpiece Adapter 25mm

NP0072



Star Handpiece Driver 25mm

RP0018

*For Angeld TMA/TMA prosthetic screw



Star Handpiece Drive 30mm

RP0019



TMA/FC Hex Handpiece Adapter 20mm

RS6194

*For Angeld TMA/TMA prosthetic screw

*For straight TMA/Flat connection




TMA/FC Hex Handpiece Adapter 25mm

RS6195A

*For straight TMA/Flat connection





UNP CloseFit™

Ultra Narrow Platform



Recommended insertion torque up to 50 Ncm

CloseFit™ UNP	Ø	L	D3	Catalog No.
	2.75	10.00	2.00	ISPUNP1027
		11.50	2.00	ISPUNP1127
		13.00	2.00	ISPUNP1327
		16.00	1.90	ISPUNP1627
		18.00	1.90	ISPUNP1827


Radiographic Template UNP CloseFit™ | ASD0615

Drill Sequence

Recommended drilling speed up to 1200 rpm

Implant Diameter

ø2.75



All measurements are in mm.

*Products may be ordered separately as an addition to the kit

Cover Screw & Healing Abutments

Recommended Torque 15(Ncm)

UNP CloseFit™ Conical Cover Screw



UNP0073

UNP CloseFit™ Slim Healing Abutment



UNP0046, GH2/L6.9

UNP0047, GH3/L7.9

UNP0048, GH4/L8.9

UNP0049, GH5/L9.9

UNP0050, GH6/L10.9

UNP CloseFit™ Healing Abutment



UNP0004, GH2/L7.1

UNP0008, GH3/L8.1

UNP0009, GH4/L9.1

UNP0012, GH5/L10.1

UNP0013, GH6/L11.1

UNP CloseFit™ Conical Healing Abutment



UNP0090, GH2/L6.9

UNP0091, GH3/L7.9

UNP0092, GH4/L8.9



UNP0093, GH5/L9.9





UNP0094, GH6/L10.9





UNP™ Impression

Recommended Torque 15(Ncm)




Transfers		L	∅	D1	H	GH	RH	Catalog No.
	UNP Slim Open Tray Transfer	18.30	5.00	2.27	14.00	-	-	UNP5614
	UNP Open Tray Transfer	18.30	6.00	2.27	14.00	-	-	UNP5514
	UNP Slim Closed Tray Transfer	15.30	3.20	2.27	13.00	-	-	UNP5113
	UNP Closed Tray Transfer	15.30	4.00	2.27	13.00	-	-	UNP5213
	UNP ScanBody	9.3	-	2.27	7.00	-	-	UNP0088
		12.3	-	2.27	10.00	-	-	UNP0089

Analog/Digital Analog		L	∅	D1	H	GH	RH	Catalog No.
	UNP Implant Analog	12.60	2.75	2.27	-	-	-	UNP0007
	UNP Digital analog 15(Ncm)	10.00	3.55	2.27	-	-	-	UNP1000



Cement Retained Abutments

Recommended Torque 25(Ncm)





Abutments	L	Ø	D1	D2	H	GH	RH	Catalog No.
	11.00	3.50	2.27	2.89	8.70	1.00	7.30	UNP0010
	9.85	4.20	2.27	3.00	7.55	1.00	6.50	UNP0019
	10.85	4.20	2.27	3.00	8.55	2.00	6.50	UNP0011
	11.85	4.20	2.27	3.00	9.55	3.00	6.50	UNP0064
	10.00	3.50	2.27	1.60	7.70	1.00	6.10	UNP0018
	9.85	4.20	2.27	2.40	7.55	1.00	6.50	UNP0063
	10.85	4.20	2.27	2.40	8.55	2.00	6.50	UNP0015
	11.85	4.20	2.27	2.40	9.55	3.00	6.50	UNP0026
	10.30	4.00	2.27	2.80	10.00	1.00	7.00	UNP0080
	11.30	4.00	2.27	2.80	11.00	2.00	7.00	UNP0081
	12.30	4.00	2.27	2.80	12.00	3.00	7.00	UNP0082
	15.10	3.50	2.27	2.85	12.80	1.10	9.30	UNP3812
	16.00	3.50	2.27	2.85	13.70	2.00	9.30	UNP3813
	17.00	3.50	2.27	2.85	14.70	3.00	9.30	UNP3814
	12.30	3.70	2.27	3.20	10.00	2.00	8.00	UNP0043
	12.30	3.70	2.27	3.20	10.00	2.00	8.00	UNP0044
	12.30	3.70	2.27	3.00	10.00	2.00	8.00	UNP0045




Overdenture Restoration





Recommended Torque 25(Ncm)



Overdentures		L	Ø	D1	D2	H	GH	RH	Catalog No.
	UNP CloseFit™ Ball Attachment	7.75	2.40	2.27	2.50	3.00	0.50	2.45	UNP0020
		8.25	2.40	2.27	2.50	3.50	1.00	2.45	UNP0021
		9.25	2.40	2.27	2.50	4.50	2.00	2.45	UNP0022
		10.25	2.40	2.27	2.50	5.50	3.00	2.45	UNP0023
		11.25	2.40	2.27	2.50	6.50	4.00	2.45	UNP0024
		12.25	2.40	2.27	2.50	7.50	5.00	2.45	UNP0025

Caps		L	Ø	D1	D2	H	GH	RH	Catalog No.
	Plastic Ball Cap White (Nylon)	2.60	4.10	-	3.60	-	-	-	RS2660
	Plastic Ball Cap Pink Soft (Nylon)	2.60	4.10	-	3.60	-	-	-	RS2662
	Titanium Ball Cap	2.90	4.70	-	4.70	-	-	-	RS2670
	Stainless Steel Ball Cap	2.90	4.70	-	4.70	-	-	-	RS2675SS

Replica		L	Ø	D1	D2	H	GH	RH	Catalog No.
	Ball Attachment Abutment Replica	16.60	4.00	2.5	-	-	0.80	2.7	RS5740

		L	Ø	Catalog No.
	UNP Abutment Screw	7.90	2.00	UNP0001
	UNP Closed Tray Transfer Screw	20.00	2.80	UNP5127
	UNP Open Tray Transfer Screw	28.00	3.10	UNP5014
	NP Retrieval Screw for Abutments	13.4	2.00	NP0014





Hand, Torque Drivers & Handpiece Drivers

Recommended Torque 25(Ncm)



Star Hand Driver -
Short 23.40L

RP0064



Star Hand Driver -
Long 28.40L

RP0004



Torque Driver for
Ball Attachment
24.4L

RS6090



Star Torque
Driver - Short 20.4L

RP0065



Star Torque
Driver-Long 25.40L

RP0005



Wrench Ratchet
Surgical Instruments
Handle 16.70L

RA0050



UNP 2mm Hex
Handpiece Adapter
Short 23.60L

UNP0068



UNP 2mm Hex
Handpiece Adapter
Medium 27.20L

UNP0071



UNP 2mm Hex
Handpiece Adapter
Long 36.00L

UNP0074



Star Handpiece
Driver - Short 25.00L

RP0018



Star Handpiece
Driver - Long 30.00L

RP0019



UNP/NP Abutment
insertion tool 25.9L

UNP1001



Founded in 2001, Adin Dental Implant Systems is a global leader in the design, development and manufacture of dental implant solutions & on track to become the world's leading independent dental implant provider.

Every business starts with a vision. Ours is simple:
We want to create the best possible dental implant solutions
that offer uncompromising quality at an affordable price.
It's our people that help our business thrive and grow.
That's why we place so much importance on building strong
personal relationships that enable us to meet our doctors'
and distributors' needs. Understanding that their success is
our success, we're focused on providing professional, high
quality and affordable solutions and exceptional service
that help them grow.

